

dynamic – dynamic range

The dynamic or dynamic range of a CCD-camera system is a widely used term to characterize the ability of a camera system to measure or image light/dark differences. In the field of photography, the dynamic range is analogous to the contrast range. However, each manufacturer defines dynamic range differently. A distinction must be made, between the “dynamic range of a CCD image sensor”, “dynamic range of an analog-to-digital-conversion” and “usable dynamic range”.

1 dynamic range of a CCD image sensor

The dynamic range (dyn_{CCD}) is defined as the ratio of the maximum possible signal (full well capacity), versus the total noise signal (in the dark). The data is expressed in decibels [dB] or is dimensionless:

$$\text{dyn}_{\text{CCD}} = \frac{\text{fullwell capacity [e}^-]}{\text{readout noise [e}^-]}$$

$$\text{dyn}_{\text{CCD}} = 20 \cdot \log\left(\frac{\text{fullwell capacity [e}^-]}{\text{readout noise [e}^-]}\right) [\text{dB}]$$

table 1: common CCD image sensor specifications

parameter	Sony ICX285AL	Kodak KAI-1020	Kodak KAI-11000
fullwell capacity [e ⁻]	18 000	40 000	60 000
readout noise [e ⁻]	5	9	11
dynamic range	3 600 : 1	4 444 : 1	5 454 : 1
dynamic range [dB]	71.1	72.9	74.7

Upon examination of these specifications, the ICX285 could be digitized using 4096 steps (corresponds to $\{2^{12}\}$ 12 bit resolution), while the KAI-11000 should be digitized using 8192 steps (corresponds to 13 bit resolution). This references to maximum values which are explained in detail in the usable dynamic range section.

Furthermore, additional gain makes sense only if the dynamic range of the A/D-converter is smaller than the dynamic range of the CCD image sensor. For example, if the dynamic range of the CCD is 70 dB, and the dynamic range of the A/D-converter is 48 dB, a useful gain or amplification would be 22 dB. Alternatively, a programmable gain would increase the performance of such a system.

2 digitization or A/D conversion dynamic range

Generated charge carriers are usually converted into voltage signals through an optimized read out circuit. These signals are amplified and finally digitized by

an analog-to-digital (A/D) converter. Therefore, the light signals (photons) are converted into digital values. The analog-to-digital converters have a given resolution or dynamic range, that in most cases is given as a power of base 2, (2^x). This means that an 8 bit resolution corresponds to 256 steps, which can be used to subdivide or convert the full scale voltage signal. CCD camera manufacturers usually optimize a combination of the dynamic range of the CCD image sensor, gain and conversion factor (defined as average conversion ratio that takes x electrons to generate one count in the image) to optimize the dynamic range of the CCD image sensor by the analog-to-digital converter.

The dynamic range of the digitization is not identical to the usable dynamic range.

table 2: dynamic range of binary resolution

resolution [bit] x => 2 ^x	dynamic range A/D conversion [digitizing steps]	dynamic range A/D conversion [dB]
8	256	48.2
10	1 024	60.2
12	4 096	72.3
14	16 384	84.3
16	65 536	96.3

The resolutions above directly correspond to the theoretical maximum limit of the converter devices. Analog-to-digital converters have an average conversion uncertainty of 0.4 - 0.7 bit, which reduces the resolution for practical applications by 1 bit. If the camera system is not limited in its dynamic range by A/D converter discrepancies, it is useful to inflate the A/D converter resolution by 1 or 2 bits. This is achieved by electronically adding a minor amount to the signal a so-called offset signal, so that the lower limit is solely provided via the image sensor and read out amplifier noise, with some resolution sacrifices from the converter.

dynamic – dynamic range

3 usable dynamic range

The actual usable dynamic range of a camera system is determined by the general settings or adjustments of the system (e.g., the quality of the dynamic range of the CCD sensor, and the compatibility between the read out circuit's dynamic range or digitization). For example, if a full well capacity of 18 000 electrons with a CCD image sensor dynamic range ratio of 4500:1 is given, getting optimal range with a 12 bit A/D converter (4096 gray levels) is possible with a conversion factor of 4.4 electrons per step (count). For very low noise cameras, it is sometimes advantageous to shift the signal resolution towards the lower noise limit and ignore the total dynamic range of the CCD image sensor, particularly if it is of minor interest in the specific application.

A particular application can also influence the usable dynamic range. Although some manufacturers attempt to define their dynamic ranges as being independent of application, these applications have a direct effect on the limits of the dynamic range. For example, in medium to high light level applications, which are photon noise limited, the intrascene dynamics may be smaller than the system's dynamic range.

4 comparing 12 bit and 8 bit dynamic range

This low light application example was recorded using a 12 bit cooled CCD camera system. The full scale, 12 bit image is completely dark. However, when the image is scaled to the 54-76 counts range, the image becomes visible. Although noise can be seen due to a weak light signal, the image can be clearly seen. Under similar conditions, but recorded with an 8 bit system

and scaling, the image can only be processed by the human eye, using the unique capabilities of the human brain, and could not be image processed in a computer.




Figure 1
Weakly illuminated image with 12bit total dynamic displayed with 255 gray levels. scaled: 0-4095 -> 0-255 graph shows readout line




Figure 2
Same image as above with 12bit total dynamic displayed with 255 gray levels. scaled: 53-70 -> 0-255 graph shows readout line




Figure 3
Same image as above but with 8bit total dynamic displayed with 255 gray levels. scaled: 7-9 -> 0-255 graph shows readout line

europe

PCO AG
Donaupark 11
93309 Kelheim, Germany

fon +49 (0)9441 2005 50
fax +49 (0)9441 2005 20
info@pco.de
www.pco.de

america

The Cooke Corporation
6930 Metroplex Drive
Romulus, Michigan 48174, USA

tel (248) 276 8820
fax (248) 276 8825
info@cookecorp.com
www.cookecorp.com