

# pc<sub>o</sub>.panda 4.2

ultra compact **sCMOS** camera

lightsheet  
scanning mode

up to 80 %  
quantum efficiency

available in  
**mono and color**

**USB 3.1**  
interface

**resolution**  
2048 x 2048 pixels  
with 6.5  $\mu\text{m}$  pixel size



65 mm

ultra  
**compact**  
design

**single cable solution**  
data & power supply via USB 3.1

# pc<sub>o</sub>.

An Excelitas Technologies Brand

## » sCMOS image sensor

<b>type of sensor</b>	scientific CMOS (sCMOS) monochrome or color (bayer pattern)
<b>resolution (h x v)</b>	2048 x 2048 active pixels
<b>pixel size (h x v)</b>	6.5 $\mu\text{m}$ x 6.5 $\mu\text{m}$
<b>sensor format / diagonal</b>	13.3 mm x 13.3 mm / 18.8 mm
<b>shutter mode</b>	rolling shutter (RS) additional feature: <b>lightsheet scanning mode<sup>1</sup></b>
<b>MTF</b>	76.9 lp/mm (theoretical)
<b>fullwell capacity</b>	45,000 e <sup>-</sup>
<b>readout noise (typ.)<sup>2</sup></b>	2.1 med e <sup>-</sup> / 2.3 rms e <sup>-</sup>
<b>dynamic range (typ.)</b>	21400 : 1 87 dB
<b>quantum efficiency</b>	up to 80 % (monochrome)
<b>spectral range</b>	370 nm ... 1100 nm
<b>dark current (typ.)</b>	15 e <sup>-</sup> /pixel/s @ 21 °C ambient temperature
<b>DSNU</b>	0.5 rms e <sup>-</sup>
<b>PRNU</b>	0.6 %
<b>anti blooming factor<sup>3</sup></b>	> 10 000

<sup>1</sup> Selectable via SDK (software development kit).

<sup>2</sup> The readout noise values are given as median (med) and root mean square (rms) values, due to the different noise models, which can be used for evaluation. All values are raw data without any filtering.

<sup>3</sup> Based on image sensor datasheet.

## lightsheet scanning mode

The PCO lightsheet scanning mode is a special readout mode dedicated to lightsheet microscopy. It is based on the rolling shutter mode in which the readout direction of the sensor is from top to bottom.


The standard line time value is 12  $\mu\text{s}$  and it can be set from this camera-specific line time up to 2 ms. Compared to the standard operation mode, the lightsheet scanning mode enables the selection of the parameters "Line Time" and "Exposure Lines". This guarantees an optimized synchronization to an existing lightsheet setup which has no selectable speed or timing. It is possible to set a delay prior to the exposure start ("delay lines").

For more information on the corresponding SDK functions, please read our pco.sdk instruction manual.

selectable parameter:  
(only via SDK)

**T<sub>L</sub> Line Time** (12  $\mu\text{s}$  ... 2 ms)  
**L<sub>E</sub> Exposure Lines** (1 ... 2048)

**Start Exposure**  
**End Exposure**



» camera system

maximum frame rate @ full resolution	40 fps
exposure / shutter time	21 $\mu$ s .. 5 s
dynamic range A/D <sup>4</sup>	16 bit
A/D conversion factor	0.65 e <sup>-</sup> /DN
pixel scan rate	44.0 MHz
pixel data rate	176.0 Mpixel/s
binning horizontal	x1, x2, x4
binning vertical	x1, x2, x4
region of interest (ROI)	horizontal: steps of 32 pixels vertical: steps of 8 pixels
non linearity	< 0.6 %
cooling method	passive cooled
trigger input signals	frame trigger, acquire (SMA connectors)
trigger output signals	exposure, busy (SMA connectors)
data interface	USB 3.1 Gen 1
time stamp	in image (1 $\mu$ s resolution)

<sup>4</sup> The high dynamic signal is simultaneously converted at high and low gain by two 12 bit A/D converters and the two 12 bit values are sophisticatedly merged into one 16 bit value.





» general

power delivery	power over USB 3.1 Gen 1
power consumption	typ. 4.5 W (max. 6.0 W)
weight	420 g
operating temperature	+ 10 °C ... + 40 °C
operating humidity range	10 % .. 80 % (non-condensing)
storage temperature range	- 10 °C .. + 60 °C
optical interface	C-mount (optional: F-mount)
maximum cable length	5 m
CE / FCC certified	yes


» frame rate table

2048 x 2048	40 fps
2048 x 1024	80 fps
2048 x 512	161 fps
2048 x 256	303 fps
2048 x 128	528 fps
1920 x 1080	76 fps
1600 x 1200	69 fps
1280 x 1024	80 fps
640 x 480	171 fps
320 x 240	321 fps

» quantum efficiency




» dimensions



F-mount and C-mount lens adapter are changeable. All dimensions are given in millimeter.


» camera view



## » applications

brightfield microscopy | fluorescence microscopy | digital pathology | single molecule localization microscopy | lightsheet fluorescence microscopy (LSFM) | calcium imaging | FRET | FRAP | structured illumination microscopy (SIM) | high-speed bright field ratio imaging | high-throughput screening | high-content screening | biochip reading | TIRF microscopy | spinning disk confocal microscopy | ophthalmology | industrial quality inspection

## » software



With pco.camware you control all camera settings, the image acquisition, and the storage of your image data. The pco.sdk is the complementary software development kit. It includes dynamic link libraries for user customization and integration on Windows PC platforms. Drivers for popular third party software packages are also available for you.

All these items like pco.camware, pco.sdk, and third party drivers are free-to-download at [www.pco.de](http://www.pco.de)

## » third party integrations



# pco.

An Excelitas Technologies Brand

telephone:

+ 49 (0) 9441 2005 50

fax:

+ 49 (0) 9441 2005 20

postal address:

Excelitas PCO GmbH

Donaupark 11

93309 Kelheim, Germany

email:

pco@excelitas.com

web:

www.pco.de

www.excelitas.com



**TII 東京インスツルメンツ**  
**TOKYO INSTRUMENTS**

グローバルにネットワークを広げ、最先端の科学をお客様に提供

本 社: 〒134-0088 東京都江戸川区西葛西6-18-14 T.I.ビル

Tel. 03-3686-4711

営業所: 〒532-0003 大阪府大阪市淀川区宮原4-1-46 新大阪北ビル

Tel. 06-6393-7411

URL: <https://www.tokyoinst.co.jp> Mail: [sales@tokyoinst.co.jp](mailto:sales@tokyoinst.co.jp)

**TII Group Company**

**UNISOKU**  
TII Group

超高真空・極低温走査型プローブ顕微鏡  
高速分光測定装置、クライオスタット

**LOTIS TII**

Nd:YAGレーザー、Ti:Sレーザー  
OPOレーザー

- 本カタログに記載されている内容は、改良のため予告無く変更する場合があります。(製品の仕様、性能、価格などはカタログ発行当時のものです)
- 本カタログに記載されている内容の一部または全部を無断で転載することは禁止されております。
- 本カタログに記載されているメーカー名、製品名などは各社の商標または登録商標です。