

data sheet  
**pco.**edge 10 bi


the next level **sCMOS** camera

**bi** back  
illuminated

resolution  
**10.4 MPixel**

pixel size  
**4.6  $\mu\text{m}$  x 4.6  $\mu\text{m}$**

interface  
**CLHS FOL**



high dynamic range  
15,385 : 1

high frame rate  
122 fps

high resolution  
4416 x 2368 pixel

low readout noise  
0.7  $e^-$ <sub>(med)</sub>

back-illuminated sCMOS  
with high MTF

temperature-stabilized  
image sensor

**pco.**

An Excelitas Technologies Brand

## technical data

### image sensor

sensor technology	back-illuminated scientific CMOS (bi sCMOS)
color type	monochrome
resolution (horizontal x vertical)	4416 pixel x 2368 pixel
pixel size (horizontal x vertical)	4.6 $\mu\text{m}$ x 4.6 $\mu\text{m}$
sensor size (horizontal x vertical)	20.3 mm x 10.8 mm
sensor diagonal	23.0 mm
shutter type	rolling shutter
modulation transfer function (theoretical max.)	108.6 lp/mm
fullwell capacity	20,000 e <sup>-</sup> @ fast scan
readout noise (typ.)	1.3 e <sup>-</sup> rms @ fast scan 1.3 e <sup>-</sup> med @ fast scan 0.8 e <sup>-</sup> rms @ slow scan 0.7 e <sup>-</sup> med @ slow scan
dynamic range (intra-scene)	15,385:1 (83.7dB) @ fast scan
peak quantum efficiency	85 % @ 500 nm
spectral range	400 nm - 1100 nm
dark current	0.4 e <sup>-</sup> /pixel/s @ +10 °C sensor temperature

### frame rate table

#### vertical resolution reduction

	fast scan	slow scan
4416 x 2368	122 fps	30 fps
4416 x 2048	141 fps	35 fps
4416 x 1024	281 fps	70 fps
4416 x 512	557 fps	139 fps
4416 x 256	1098 fps	274 fps
4416 x 128	2132 fps	533 fps
4416 x 64	4028 fps	1007 fps
4416 x 32	7252 fps	1813 fps
4416 x 16	12,086 fps	3021 fps
4416 x 8	18,130 fps	4532 fps

#### typical resolutions


	fast scan	slow scan
2304 x 2304	125 fps	31 fps
2048 x 2048	141 fps	35 fps
1920 x 1080	266 fps	66 fps
1280 x 1024	281 fps	70 fps
640 x 512	557 fps	139 fps
320 x 256	1098 fps	274 fps

## camera

<b>max. frame rate @ full resolution</b>	122 fps @ fast scan 30 fps @ slow scan
<b>exposure time range</b>	6.8 $\mu$ s - 1 s @ fast scan 27.5 $\mu$ s - 1 s @ slow scan
<b>dynamic range A/D</b>	16 bit
<b>conversion factor<sup>1</sup></b>	0.275 e <sup>-</sup> /DN @ fast scan 0.305 e <sup>-</sup> /DN @ slow scan
<b>pixel rate</b>	1467 MPixel/s @ fast scan 366 MPixel/s @ slow scan
<b>region of interest (ROI)</b>	horizontal: steps of 1 column vertical: steps of 8 rows
<b>binning</b>	horizontal: x2, x4 (average, sum) vertical: x2, x4 (average, sum)
<b>non-linearity</b>	< 0.33 % @ fast scan < 0.5 % @ slow scan
<b>dark signal non-uniformity (DSNU)</b>	< 0.23 e <sup>-</sup> rms @ fast scan < 0.07 e <sup>-</sup> rms @ slow scan
<b>photo response non-uniformity (PRNU)</b>	< 0.3 %
<b>cooling temperature image sensor</b>	+10 °C stabilized
<b>cooling method</b>	forced air & water
<b>trigger input signals</b>	frame trigger, sequence trigger, programmable input
<b>trigger output signals</b>	exposure, busy, programmable output
<b>input / output signal interface</b>	SMA connectors
<b>time stamp</b>	in image (1 $\mu$ s resolution)
<b>data interface</b>	Camera Link HS FOL

<sup>1</sup> According to EMVA1288 the conversion factor equals the inverse of the system gain and can be operational mode dependent.

## quantum efficiency



## general


power supply	24 VDC ( $\pm 10\%$ )
power consumption	< 40 W
weight	1.35 kg
dimensions (height x width x length)	95 mm x 90 mm x 109 mm
operating temperature range	+10 °C to +40 °C
operating humidity range (non-condensing)	10 % to 80 % (non-condensing)
storage temperature range	-10 °C to +60 °C
CE / FCC certified	yes

## optical interface

direct mounting	6.2 mm $\pm 10\%$
lens mounting	C-Mount
optional lens mounting	F-Mount, TFL-Mount
optional lens remote controller	EF-Mount, EF-S-Mount (Canon)

Configure your optical setup with our **MachVis Lens Selector** online tool.

## dimensions



Outlines of pco.edge 10 bi (all dimensions given in mm).

## software

Our main camera control software pco.camware is the first choice to get started with your camera. It enables full control of all camera settings and makes image acquisition and storage very easy. Using different layouts, stiles and features you can customize it exactly to your needs.



### You are using a different software:

PCO cameras are also integrated in a variety of software applications. Check our homepage to find a list of all applications that support PCO cameras.






(in preparation)



(in preparation)

### You want to create your own application for the camera:

We offer a wide range of Software Development Kits (SDK) for different programming languages, both for windows and linux. Our pco.sdk, pco.recorder and high-level SDK are designed for C/C++ apps. With pco.python, pco.matlab, pco.labview and pco.java you can control the camera in your C#, python, matlab, labview and java applications, respectively.



### Your use case is in the field of microscopy:

PCO cameras are also integrated in µManager.



## areas of application

bright-field microscopy | fluorescence microscopy | digital pathology | single molecule localization microscopy (SMLM) | lightsheet fluorescence microscopy (LSFM) | selective plane illumination microscopy (SPIM) | structured illumination microscopy (SIM) | raman spectroscopy | calcium imaging | Förster resonance energy transfer (FRET) | fluorescence recovery after photobleaching (FRAP) | high-speed bright-field ratio imaging | high throughput screening | ophthalmology | biochip reading | total internal reflection fluorescence microscopy (TIRF) | 3D metrology | industrial quality inspection | wafer inspection | image intensifier imaging | intravital microscopy | inspection | material testing | biometrics | in-vivo microscopy

# pco.

An Excelitas Technologies Brand

postal address:	Excelitas PCO GmbH Donaupark 11 93309 Kelheim, Germany
telephone:	+49 (0) 9441 2005 0
e-mail:	<a href="mailto:pco@excelitas.com">pco@excelitas.com</a>
web:	<a href="http://www.excelitas.com/pco">www.excelitas.com/pco</a>



**EXCELITAS**  
TECHNOLOGIES®



グローバルにネットワークを広げ、最先端の科学をお客様に提供

本社：〒134-0088 東京都江戸川区西葛西6-18-14 T.I.ビル

Tel.03-3686-4711

営業所：〒532-0003 大阪府大阪市淀川区宮原4-1-46 新大阪北ビル

Tel.06-6393-7411

URL：<https://www.tokyoinst.co.jp> Mail：sales@tokyoinst.co.jp



超高真空・極低温走査型プローブ顕微鏡  
高速分光測定装置、クライオスタット



Nd:YAGレーザー、Ti:SLレーザー  
OPOLレーザー

- 本カタログに記載されている内容は、改良のため予告無く変更する場合があります。（製品の仕様、性能、価格などはカタログ発行当時のものです）
- 本カタログに記載されている内容の一部または全部を無断で転載することは禁止されております。
- 本カタログに記載されているメーカー名、製品名などは各社の商標または登録商標です。