

Description

The PhotoniQ model MCPC618 is a complete, off-the-shelf, high speed, eight channel pulse counting system for PMTs, silicon photomultipliers (SiPM) and APDs. Implemented as a stand-alone laboratory instrument with a PC interface, the MCPC618 is used for discrimination, counting, and data acquisition (DAQ) of single pulse events across eight independent counting channels. Its unique front end design permits direct connection to most PMTs without the need for an external preamplifier. Flexible intelligent triggering allows the unit to reliably acquire count data using one of several sophisticated triggering techniques. The MCPC618 is fully configurable through the PC via its USB 2.0 port using an included graphical user interface. Continuous high speed data transfers to the PC are also handled through this port. Additionally, a LabView™ generated DLL is provided for users who wish to write their own applications that interface directly to the unit.


Applications

- Fluorescence Spectroscopy
- Fluorescence Lifetime Measurement
- Chemiluminescence Detection
- Bioluminescence Detection
- Photon Correlation Spectroscopy
- DNA Sequencing
- Particle Sizing
- LIDAR
- Flow Cytometry
- Neutrino Detection
- Spatial Radiation Detection
- Confocal Microscopy
- Particle Physics



Features

- Includes eight independent AC-coupled counting channels with on-board discriminators
- Maximum count rate greater than 250 MHz per channel for a total of 2 billion counts per second
- Pulse pair resolution of less than 4 nsec
- Adjustable discriminator threshold control
- Flexible control of counting period parameters such as delay, width, or external boxcar
- Adjustable microGate provides additional level of gating at sub-nanosecond time resolution
- Synchronization of microGate to external excitation source
- Parallel, high speed hardware processor unit performs real-time data filtering and dark count subtraction
- Programmable data filtering function for real time detection of predefined energy patterns or spectrums
- Advanced triggering supports external, internal, level, boxcar, pre-triggering and input threshold crossing modes
- Trigger stamping and time stamping with 100 nsec resolution
- USB 2.0 interface supports high data transfer rates
- Graphical User Interface (GUI) for menu driven data acquisition and configuration
- LabVIEW™ generated DLL for interface to custom applications



Software Features & Functions

- Graphical User Interface (GUI) for menu driven data acquisition, configuration, and status
- Real time display shows total count for count period across all channels for each trigger
- Integrated log file viewer permits on-screen viewing of logged count records
- Microgate width and delay controllable in GUI with 500 psec resolution
- Fully programmable count period minimizes dead time
- High speed trigger counter and record counter
- Acquisition can be programmed to acquire for a preset number of records
- Trigger stamping and record time stamping with 100 nsec resolution
- USB 2.0 interface supports high transfer rates
- Included Microsoft Windows DLL for interface to custom user applications


Included Accessories and Software

The MCPC618 is enclosed in a rugged, EMI-shielded, instrument case and is shipped with the following standard components and software:

- PhotoniQ Control and Acquisition Interface Software CD-ROM
- DC power supply (+5V, 3A) with power cord
- USB 2.0 cable

Specifications

Recommended Sensor Interface Circuits



Description	Specification
Number of Channels	8
Input Impedance	50 ohm
AC-Coupling Time Constant	5 usec
Pulse Pair Resolution (PPR)	4 nsec max.
Minimum Detectable Pulse Amplitude	8 mV
Maximum Count Rate per Channel	250 MHz
Count Period Range	100 nsec to 1 sec
Maximum Count per Count Period	16,383
Maximum Trigger Rate	400 KHz
Sustained Trigger Rate (8 Channels Enabled)	250 KHz
Power Consumption	3.5 Watts typ.
Width	9.843 in. (250 mm)
Height	3.346 in. (85 mm)
Length	10.236 in. (260 mm)

* See MCPC618 User Manual for details



Vertilon Corporation has made every attempt to ensure that the information in this document is accurate and complete. Vertilon assumes no liability for errors or for any incidental, consequential, indirect, or special damages including, without limitation, loss of use, loss or alteration of data, delays, lost profits or savings, arising from the use of this document or the product which it accompanies.

Vertilon reserves the right to change this product without prior notice. No responsibility is assumed by Vertilon for any infringements of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under the patent and proprietary information rights of Vertilon Corporation.

© 2019 Vertilon Corporation, ALL RIGHTS RESERVED

PS2720.1.2. Jul 2019

Vertilon Corporation, 66 Tadmuck Road, Westford, MA 01886 / Tel: (978) 692-7070 / Fax: (978) 692-7010 / www.vertilon.com

Description

The PhotoniQ model MCPC680 is a complete, off-the-shelf, high speed, 32 channel pulse counting system for PMTs, silicon photomultipliers (SiPM) and APDs. Implemented as a stand-alone laboratory instrument with a PC interface, the MCPC680 is used for discrimination, counting, and data acquisition (DAQ) of single pulse events across 32 independent counting channels. Its front end design permits direct connection to most PMTs and SiPMs without the need for a preamplifier. Flexible intelligent triggering allows the unit to reliably acquire count data using one of several sophisticated triggering techniques. The MCPC680 is fully configurable through the PC via its USB 2.0 port using an included graphical user interface. Continuous high speed data transfers to the PC are also handled through this port. Additionally, a LabVIEW™ generated DLL is provided for users who wish to write their own applications that interface directly to the unit.


Applications

- Fluorescence Spectroscopy
- Fluorescence Lifetime Measurement
- Chemiluminescence Detection
- Bioluminescence Detection
- Photon Correlation Spectroscopy
- Bioaerosol Detection and Discrimination
- LIDAR
- Optical Tomography of Biological Tissue
- Low Light Level Detection
- Flow Cytometry
- Single Molecule Detection
- Neutrino Detection
- Spatial Radiation Detection
- Confocal Microscopy
- Particle Physics



Features

- Includes 32 independent counting channels with PMT / SiPM interface and discriminators
- Threshold control of internal leading / trailing edge discriminators
- Pulse width discrimination based on user-defined minimum / maximum conditions
- Pulse pair resolution of less than 4 nsec.
- Maximum count rate greater than 250 MHz per channel for a total of 8 billion counts per second
- Intelligent triggering supports external, internal, level, and continuous counting modes
- Flexible control of counting period parameters such as delay, width, or external
- Adjustable microGate provides additional level of count gating at sub-nanosecond time resolution
- Synchronization of microGate to external excitation source
- Trigger stamping and time stamping with 100 nsec resolution
- USB 2.0 interface supports high data transfer rates
- Graphical User Interface (GUI) for menu driven data acquisition and configuration
- LabVIEW™ generated DLL for interface to user custom applications



Software Features & Functions


- Graphical User Interface (GUI) for menu driven data acquisition, configuration, and status
- Real time display shows total count for count period across all channels for each trigger
- Integrated log file viewer permits on-screen viewing of logged count records
- Microgate width and delay controllable in GUI with 500 psec resolution
- Fully programmable count period minimizes dead time
- High speed trigger counter and record counter
- Acquisition can be programmed to acquire for a preset number of records
- Trigger stamping and record time stamping with 100 nsec resolution
- USB 2.0 interface supports high transfer rates
- Included Microsoft Windows DLL for interface to custom user applications

Included Accessories and Software

The MCPC680 is enclosed in a rugged, EMI-shielded, instrument case and is shipped with the following standard components and software:

- PhotonIQ Control and Acquisition Interface Software CD-ROM
- DC power supply (+5V, 3A) with power cord
- USB 2.0 cable

Recommended Sensor Interface Circuits



Specifications

Description	Specification
Number of Channels	32
Input Impedance	50 ohm, AC coupled
Pulse Pair Resolution (PPR)	4 nsec max.
Minimum Detectable Pulse Amplitude	8 mV
Maximum Count Rate per Channel	250 MHz
Count Period Range	50 nsec to 1 sec
Maximum Count per Count Period	16,383
Maximum Trigger Rate	125 KHz
Sustained Trigger Rate (32 Channels Enabled)	64 KHz
Power Consumption	5 Watts typ.
Width	9.843 in. (250 mm)
Height	3.346 in. (85 mm)
Length	10.236 in. (260 mm)

* See MCPC680 User Manual for details



Vertilon Corporation has made every attempt to ensure that the information in this document is accurate and complete. Vertilon assumes no liability for errors or for any incidental, consequential, indirect, or special damages including, without limitation, loss of use, loss or alteration of data, delays, lost profits or savings, arising from the use of this document or the product which it accompanies.

Vertilon reserves the right to change this product without prior notice. No responsibility is assumed by Vertilon for any infringements of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under the patent and proprietary information rights of Vertilon Corporation.

© 2020 Vertilon Corporation, ALL RIGHTS RESERVED

PS2741.1.1 Aug 2020

Description

The PhotoniQ model MCPC682 is a complete, off-the-shelf, high speed, 64 channel pulse counting system for PMTs, silicon photomultipliers (SiPM) and APDs. Implemented as a stand-alone laboratory instrument with a PC interface, the MCPC682 is used for discrimination, counting, and data acquisition (DAQ) of single pulse events across 64 independent counting channels. Its front end design permits direct connection to most PMTs and SiPMs without the need for a preamplifier. Flexible intelligent triggering allows the unit to reliably acquire count data using one of several sophisticated triggering techniques. The MCPC682 is fully configurable through the PC via its USB 2.0 port using an included graphical user interface. Continuous high speed data transfers to the PC are also handled through this port. Additionally, a LabVIEW™ generated DLL is provided for users who wish to write their own applications that interface directly to the unit.


Applications

- Fluorescence Spectroscopy
- Fluorescence Lifetime Measurement
- Chemiluminescence Detection
- Bioluminescence Detection
- Photon Correlation Spectroscopy
- Bioaerosol Detection and Discrimination
- LIDAR
- Optical Tomography of Biological Tissue
- Low Light Level Detection
- Flow Cytometry
- Single Molecule Detection
- Neutrino Detection
- Spatial Radiation Detection
- Confocal Microscopy
- Particle Physics




Features

- Includes 64 independent counting channels with PMT / SiPM interface and discriminators
- Threshold control of internal leading / trailing edge discriminators
- Pulse width discrimination based on user-defined minimum / maximum conditions
- Pulse pair resolution of less than 4 nsec.
- Maximum count rate greater than 250 MHz per channel for a total of 16 billion counts per second
- Intelligent triggering supports external, internal, level, and continuous counting modes
- Flexible control of counting period parameters such as delay, width, or external
- Adjustable microGate provides additional level of count gating at sub-nanosecond time resolution
- Synchronization of microGate to external excitation source
- Trigger stamping and time stamping with 100 nsec resolution
- USB 2.0 interface supports high data transfer rates
- Graphical User Interface (GUI) for menu driven data acquisition and configuration
- LabVIEW™ generated DLL for interface to user custom applications



Software Features & Functions


- Graphical User Interface (GUI) for menu driven data acquisition, configuration, and status
- Real time display shows total count for count period across all channels for each trigger
- Integrated log file viewer permits on-screen viewing of logged count records
- Microgate width and delay controllable in GUI with 500 psec resolution
- Fully programmable count period minimizes dead time
- High speed trigger counter and record counter
- Acquisition can be programmed to acquire for a preset number of records
- Trigger stamping and record time stamping with 100 nsec resolution
- USB 2.0 interface supports high transfer rates
- Included Microsoft Windows DLL for interface to custom user applications

Included Accessories and Software

The MCPC682 is enclosed in a rugged, EMI-shielded, instrument case and is shipped with the following standard components and software:

- PhotonIQ Control and Acquisition Interface Software CD-ROM
- DC power supply (+5V, 3A) with power cord
- USB 2.0 cable

Recommended Sensor Interface Circuits



Specifications

Description	Specification
Number of Channels	64
Input Impedance	50 ohm, AC coupled
Pulse Pair Resolution (PPR)	4 nsec max.
Minimum Detectable Pulse Amplitude	8 mV
Maximum Count Rate per Channel	250 MHz
Count Period Range	50 nsec to 1 sec
Maximum Count per Count Period	16,383
Maximum Trigger Rate	125 KHz
Sustained Trigger Rate (64 Channels Enabled)	32 KHz
Power Consumption	5 Watts typ.
Width	9.843 in. (250 mm)
Height	3.346 in. (85 mm)
Length	10.236 in. (260 mm)

* See MCPC682 User Manual for details



Vertilon Corporation has made every attempt to ensure that the information in this document is accurate and complete. Vertilon assumes no liability for errors or for any incidental, consequential, indirect, or special damages including, without limitation, loss of use, loss or alteration of data, delays, lost profits or savings, arising from the use of this document or the product which it accompanies.

Vertilon reserves the right to change this product without prior notice. No responsibility is assumed by Vertilon for any infringements of patents or other rights of third parties which may result from its use. No license is granted by implication or otherwise under the patent and proprietary information rights of Vertilon Corporation.

© 2020 Vertilon Corporation, ALL RIGHTS RESERVED

PS2742.1.2 Aug 2020