


The modular Software platform



Preparation of your experiment

- ☞ **EductX:** Automated integration of reactant specific data from an ELN (data base).
- ☞ **OperX:** Create your recipe wherever from you want.
- ☞ **SamX:** Your assistant to switch peripheral devices easily.


Compliance, comfort and safety

- ☞ **SecureX:** Makes your software GMP compliant (CFR 21 part 11).
- ☞ **RemoteX:** Control your reactor from wherever you want.
- ☞ **MessageX:** Your email alert.

Data Management

- ☞ **ChartX:** Graphical view (trend) of your experiment.
- ☞ **e-Journal:** Automatic generated lab journal of your experiment.
- ☞ **CollectX:** Automatic data transfer to an ELN, LIMS, cloud etc.


FlexySys - The modular software platform for your laboratory application



FlexySys — simplicity and flexibility through structured functions


Simplicity: Thanks to intuitive functions, experiments can be carried out safely and without extensive training.

Flexibility: Thanks to a wide range of standardized functions, we can offer you a solution tailored to your own process, so that you can conduct your work as efficiently as possible. Existing equipment can also be integrated into the software. This way, you not only save money, but also increase the system's availability.



Efficiency, safety and reproducibility thanks to recipe control


Using recipes, sub-processes such as inerting or even complete experiments can be carried out reproducibly and efficiently, even without supervision. Maximum flexibility is guaranteed by the combination of manual interventions, fully automatic recipe operation and the “edit on the fly” function. Alongside all the necessary safety limits, which the system immediately regulates into the previously defined safety state, a variety of process limits can also be defined. These include, for example, the maximum permissible temperature rise during dosing.



Customer-specific adjustments

The software can be tailored to a large number of different processes. For example, distillations, filtrations or pressure controls can be automated via the software using standardized functions, while reaction energies can be measured (calorimetry) or analysis devices such as turbidity measurements and particle size analyzers can be implemented.

Customer-specific turnkey solutions, combined with services in the field of plant design and plant qualification in the GMP environment (IQ/OQ), protect your investments, thanks to the modular way that they can be adapted.



Data management and eJournal

During an experiment, all the events and data are recorded automatically. This also applies for any integrated analytical instrument. In addition, all the data along the workflow, such as the numbers of manual weighings of solids or the batch numbers of educts, can be managed via the software. All the data and information is compiled in Word format in an automatically generated e-journal, which can then be centrally archived in higher-level data management programs (ELN or LIMS) using the “CollectX” add-on. This way, the traceability of all the experiment related data is ensured and data analysis is also guaranteed across departments.

SYSTAG, System Technik AG

Bahnhofstr. 76 | CH-8803 Rüschiikon
 Phone +41 44 704 54 54 | Fax +41 44 704 54 55
 infos@systag.ch | www.systag.ch

TII 東京インスツルメンツ 本社: 〒134-0088 東京都江戸川区西葛西6-18-14 T.ビル Tel.03-3686-4711
TOKYO INSTRUMENTS 営業所: 〒532-0003 大阪府大阪市淀川区宮原4-1-46 新大阪北ビル Tel.06-6393-7411
 グローバルにネットワークを広げ、最先端の科学をお客様に提供 URL: https://www.tokyoinst.co.jp Mail: sales@tokyoinst.co.jp

TII Group Company **UNISOKU** 超高真空・極低温走査型プローブ顕微鏡 超高速分光測定装置、クライオスタット **LOTIS TII** Nd:YAGレーザー、Ti:Sレーザー OPOレーザー

※本カタログに記載されている内容は、改良のため予告無く変更する場合があります。(製品の仕様、性能、価格などはカタログ発行当時のものです)
 ※本カタログに記載されている内容の一部または全部を無断で転載することは禁止されております。
 ※本カタログに記載されているメーカー名、製品名などは各社の商標または登録商標です。