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
The screenshot displays the Nanosculptor TII software interface. The main window shows a 3D topographic map of a sample with a highlighted region. The interface includes a menu bar (File, View, Tools, Help), a toolbar with 3D navigation tools (XY, XZ, YZ), and a central camera view labeled '顕微鏡像' (Microscope Image). On the left, there is a 'Camera' panel with a table of shapes and their corresponding icons, and a 'Batch' panel with buttons for 'Batch', 'Tiling', 'Processing Order', 'Overlap', 'Return', 'Data Deletion', and 'Data Editing'. On the right, there is a 'Settings' panel with checkboxes for 'Stage Current Position', 'Coordinate Display', 'Stage Movable Range', 'Carrier', 'Sample Shape', 'Processing Data', 'Grid', 'Processing Area', and 'Area Coordinate Display'. Below the main view, there is a 'Stage Control' panel with directional buttons and a table of current and target values for X, Y, Z, and A. The bottom panel contains 'Processing Parameters' and 'Presets' sections, with buttons for 'Move', 'Acquire', and 'AF'. The status bar at the bottom indicates 'ステージ停止中' (Stage stopped).

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

The Notepad window displays a 3D coordinate system with X, Y, and Z axes. A starting point (始点) is marked at the origin. Three planes are defined at different Z heights: 平面① at $z=10000\text{nm}$, 平面② at $z=5000\text{nm}$, and 平面③ at $z=0\text{nm}$. Red arrows indicate the direction of processing on each plane. The text data in the Notepad window is as follows:

```
0 0 10000
0 50000 10000
50000 50000 10000
50000 0 10000
0 0 10000

0 0 5000
0 50000 5000
50000 50000 5000
50000 0 5000
0 0 5000

0 0 0
0 50000 0
50000 50000 0
50000 0 0
0 0 0
```

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The 'Create scan array' dialog box has the following settings:

Quantity of points	Value of one step	Inverse direction
X: 25	2000 nm	<input type="checkbox"/>
Y: 25	2000 nm	<input type="checkbox"/>
Z: 5	1000 nm	<input checked="" type="checkbox"/>


Buttons: Create, Cancel

Options: XYZ (dropdown), Order of points creating, Creation of lines, Optimum scanning

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The screenshot shows the NMG Data Processing software interface. The main workspace displays a grid with a blue outline of a flower-like shape on the left and a red bitmapped version of the same shape on the right. The text '複雑な形状データ' (Complex shape data) is positioned above the red shape. At the bottom of the grid, the text 'さくら2012' (Sakura 2012) is displayed, with '数値/文字データ' (Numerical/Text data) below it. The software's toolbar and control panels are visible at the top, including options for file operations, display settings, and data processing parameters.

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The screenshot displays the Nanosculptor TII software interface (Ver. 0.01) with the following components:

- Top Panel:** Menu bar (File, View, Tools, Help) and a toolbar with icons for file operations, 3D/XY/XZ/YZ views, camera, and execution.
- XY View (Top Left):** 2D grid view with a grid size of 946.25um.
- XZ View (Bottom Left):** 2D grid view with a grid size of 946.25um.
- 3D View (Bottom Right):** 3D perspective view of the grid-based workpiece with a grid size of 946.25um.
- Microscope Image (Top Right):** A live microscope image of the workpiece, labeled "顕微鏡像".
- Right Panel:** A settings sidebar with sections for Stage (Stage position, coordinate display, range), Sample (Shape), Processing Data (Data, Grid, Area, Area coordinate display), Registration Point (Reference point, origin, loading position), Color specification (Color bar), Line width, and Dot size.
- Bottom Panel:** Control and parameter sections.
 - Stage Control:** Includes directional movement buttons, a table for current, target, and step values for X, Y, Z, and A axes, and a movement speed of 10000.
 - Processing Parameters:** Includes a table for current, target, and step values for X, Y, Z, and A axes, and a "移動" (Move) button.
 - Laser Parameters:** Includes fields for Laser wavelength (800nm), Frequency (1000 Hz), Objective lens (x5), Energy (30 mW), and Laser correction, with "セット" (Set) and "再取得" (Refresh) buttons.
 - Additional Controls:** Includes buttons for "詳細設定" (Detailed Settings), "周波数設定" (Frequency Setting), "ステージドライブ" (Stage Drive), "ステージドライブ設定" (Stage Drive Setting), "パワーメータ" (Power Meter), and "AFディテクト" (AF Detect).

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