Arcoptix SA CH-2000 Neuchâtel Switzerland Phone : +41 32 731 04 66 Email: info@arcoptix.com



ARCOPTIX OEM010



The ARCoptix OEM010 is our most compact FT-IR complete system. It features a built-in SiC globar heated to 1550 °K and regulated in temperature to ensure stable illumination. The light that is being output by the instrument is first modulated by our permanently aligned interferometer, making it robust against stray light that might eventually reach the detector. The OEM010 rugged casing also holds a highly sensitive TE-MCT detector which maximizes the response of the system.



Optical routing in the ARCOPTIX OEM010

Features Internal light source Compact, rugged unit Custom I/O coupling modules on request Robust to ambient light perturbation Dynamically adjustable resolution: 8cm⁻¹ 4cm⁻¹ 2cm⁻¹ 0.5 cm⁻¹ (on request) Wear free moving parts for extended lifetime No purging of the interferometer required **Temperature controlled** reference laser

Arcoptix SA CH-2000 Neuchâtel Switzerland Phone : +41 32 731 04 66 Email: info@arcoptix.com



Specifications

Product code	FTIR-OEM010-060-4TE	FTIR-OEM010-085-4TE	FTIR-OEM010-120-4TE
Beam-splitter material	CaF ₂		ZnSe
Spectral Range [cm ⁻¹]	1'660-5'000	1'200-6'600	830-5'000
Spectral Range [µm]	2-6	1.5-8.5	2-12
Detector Type	МСТ	МСТ	
	(2-TE cooled)	(4-TE cooled)	
Detector Peak D* [cm Hz ^{1/2} W ⁻¹]	>1x10 ¹¹	>8x10 ⁹	>4x10 ⁹
Signal-to-noise ratio (SNR)	> 80'000:1 ⁱ	>40'000:1 ⁱ	
Output beam characteristics	Ø 12.7mm collimated (max ~30mrad half angle)		
Interferometer type	Permanently aligned, double retro-reflector design		
Resolution (unapodized) [cm ⁻¹]	0.5 ⁱⁱ , 2, 4, 8 (user selectable)		
Wavenumber repeatability	<10 PPM		
Scan frequency	>4 hz @ 4cm ⁻¹		
Internal reference laser	Temperature-stabilized solid-state laser @850nm		
A/D Converter	24 bit		
Amplifier	4 gain levels low noise trans-impedance amplifier		
Operating temperature [°C]	10 to 40		
Storage temperature [°C]	-20 to 60		
Power requirement	12V / 40W		
Communication Interface	USB 2.0		
Software Interface	Windows 7/10/11, Linux		
	API for controlling the instrument via our DLL		
	Single board computer available on request		
Dimensions [mm]	165x145x82		
Weight [g]	2100		

ⁱ Input to output direct light routing, 60s measurement, around peak sensitivity wavelength, 4cm-1, Norton-Beer weak apodization.

ⁱⁱ Available on request only, please contact us at <u>info@arcoptix.com</u> for details SPECIFICATIONS ARE SUBJECT TO CHANGES WITHOUT NOTICE.



グローバルにネットワークを広げ、最先端の科学をお客様に提供 URL: https://www.tokyoinst.co.jp Mail: sales@tokyoinst.co.jp

TI Group Company

