

PHOTON RT 0420 Ultra

**VIS-NIR TRANSMISSION MEASUREMENT OF CUTTING-EDGE COATINGS
WITH STEEP SLOPES AND HIGH OPTICAL DENSITY**



UNIQUE MEASUREMENT CAPABILITIES:

- Broadband measurement with one click
- Optical density: up to OD8
- Filter bandwidth: down to 0,1 nm FWHM
- Slope steepness: down to 0,02 %

SPECIFICATIONS

PARAMETER	DESCRIPTION
MODEL	0420 Ultra
OPTICAL CONFIGURATION	
Optical scheme of monochromator	Double Additive Czerny-Turner
Optics	Mirror, protected Ag
Reference channel	Yes
Wavelength sampling pitch, nm	0,01 to 100
Spot size on the measured sample, mm	4 x 3
Photometric functions	%T, OD
Effective wavelength range, nm	340 - 2000
Ultimate spectral resolution, nm 340-990 nm 990-2000 nm	0,06 0,13
Wavelength accuracy, nm 340-990 nm 990-2000 nm	+/- 0,02 +/- 0,04
Wavelength repeat accuracy, nm 340-990 nm 990-2000 nm	+/- 0,01 +/- 0,02
Ultimate optical density, -log ₁₀ (T) 340-420 nm 420-925 nm 925-2000 nm	6 8 7
Stray light level, % at 532 nm	< 0,000001
Angle of beam divergence, deg	+/- 0,2
Photometric accuracy	(VIS) NIST SRM 930e: +/-0.003 Abs (1 Abs) NIST SRM 1930: +/-0.003 Abs (0.33 Abs); +/-0.006 Abs (2 Abs)
Photometric repeat accuracy	(VIS) NIST SRM 930e: 0.0004 Abs (1 Abs) NIST SRM 1930: 0.0001 Abs (0.33 Abs); 0.005 Abs (2 Abs) Determined using 0,1 second accumulation, maximum deviation for 10 subsequent measurements
Stability of baseline, %/hour	< 0,2 (one hour warm-up time)
Light sources	High energy plasma source, HgAr wavelength calibration lamp. All sources built-in and pre-configured
SAMPLE COMPARTMENT	
Standard sample stage (included)	Measurement of transmission of plano samples with size greater than 8 x 10mm
Size of samples	Min. 6x4 mm Max. sample size: - up to 55 (90) mm x 120 mm with closed lid for standard sample stage - up to 152,4 mm (round) with closed lid for optional motorized sample stage
INTERFACE, DIMENSIONS AND WEIGHT	
Interface	USB 2.0, Windows-based, English
File saving options	res (txt), xls, pdf, csv
Power consumption, Watt	250
Power input	110-220 V, 50-60 Hz
Width x Depth x Height, mm (inches)	820 x 545 x 360 (32 2/5" x 25 2/5" x 14 1/5")
Net weight, kg (lbs)	55 (121)
Values measured after 60 minutes warm-up time	

- 本カタログに記載されている内容は、改良のため予告無く変更する場合があります。（製品の仕様、性能、価格などはカタログ発行当時のものです）
- 本カタログに記載されている内容の一部または全部を無断で転載することは禁止されております。
- 本カタログに記載されているメーカー名、製品名などは各社の商標または登録商標です。