

Solo 698.4 QT

Ultra-narrow Linewidth Laser



UniKLasers supplies ultra-narrow linewidth lasers from our QT series to our partners in quantum sensing, metrology and other technologies at the specific wavelengths related to the exact atomic transitions they target. Our QT laser series includes the Solo 813.42 QT and Solo 780.24 QT for Rubidium, Solo 698.4 QT and Solo 689.4 QT for Strontium.

KEY FEATURES



Ultra-narrow Linewidth
 ≤ 0.3 MHz



High Power Stability
 ≤ 2.0 % over 8 hours



High Spectral Stability
 ± 1 pm over 8 hours



Excellent Beam Quality



Designed for Strontium Atoms

APPLICATIONS

Quantum Research - Strontium

Specifications

Output Beam Parameters:

Output Power	10, 30, or 50 mW*
Wavelength	698.4 nm
Spectral Bandwidth	≤ 0.3 MHz
Spatial Mode	TEM ₀₀
Spectral Stability	± 1.0 pm (over 8 hour operation)
Coherence Length	> 100 m
Output Power Stability	≤ 2.0 % (over 8 hour operation)
Output Power Noise	≤ 0.1 % RMS (10 Hz – 10 MHz)
Beam Divergence	≤ 1 mrad, diffraction limited
Beam Diameter at Output Aperture	0.8 – 1.2 mm
Beam Pointing Stability	≤ 5 μ rad/°C
Polarisation Ratio	$\geq 100:1$, vertical

Integration Features:

Plug-In USB Connectivity	Combined Heatsink
Versatile Control Software	Remote Diagnostic Support

Laser Head Dimensions:

L x W	207 x 80 mm
Beam Height	75 mm

Environmental Conditions:

Ambient Temperature Range	18 – 30 °C
Laser Head Interface Stability	± 1.5 °C
Storage	0 – 50 °C
Humidity	5 – 95 %, non-condensing

Optional Accessories:


Heatsink	Fan-Assisted
	Water-Cooled with Thermoelectric Chiller
External Power Control Module	0 – 100 %, continuous
Tuning Range, mode-hop free	25 – 50 GHz

Low Power Alignment Beam Mode

* Other output powers available on request

Laser and Heatsink Dimensions

238 x 80 x 95 mm



CONTACT US

- ✉ info@uniklasers.com
- ☎ +44 (0)131 333 2200
- 🌐 uniklasers.com



<https://www.tokyoinst.co.jp/>

本社：〒134-0088 東京都江戸川区西葛西 6-18-14 T. I. ビル
TEL : 03-3686-4711 FAX : 03-3686-0831
大阪営業所：〒532-0003 大阪市淀川区宮原 4-1-46 新大阪北ビル
TEL : 06-6393-7411 FAX : 06-6393-7055

DESIGNED &
MANUFACTURED
IN THE UK



AVOID EYE OR SKIN EXPOSURE TO
DIRECT OR SCATTERED RADIATION
CLASS 4 LASER PRODUCT
MAX. POWER: 1000 mW
(IEC 60825-1)