


Pen-Ray® Accessories



Lamp Shields & Filters		
Part Number	Description	Special
98-0008-01	Lamp Shield A	0.04" (1.02mm) ID
98-0008-02	Lamp Shield B	0.31" x 0.63" (7.78 x 16mm) window
98-0008-03	Lamp Shield C	0.19" x 1.5" (4.8 x 38.1mm) window
98-0017-01	G-275 Filter, SW	Absorbs visible light while transmitting 254nm
98-0018-01	G-278 Filter, SW	Absorbs shortwave while transmitting longwave

Shields and Filters

Pen-Ray Filter and Shields, for directing or modifying the emission output, can be ordered for the 11SC-1 Lamp.

- Shield A has a 0.04 inch ID hole for a point-like source
- Shield B a 0.31 x 0.63" window;
- Shield C a 0.19 x 1.5" window
- G-275 filter absorbs visible light while transmitting ultraviolet at 254nm
- G-278 filter converts shortwave radiation to longwave at 365nm