


The Most Advanced Vibration Absorption Systems

SOTOMOXXX - SOTO MODULAR

A lightweight and compact vibration absorption systems allows to significantly reduce unwanted vibrations translated to working environment from various environment disturbances. Developed by engineers, physicists, and leading industry consultants with expertise in vibration isolation and precision instruments, SDTDMD effectiveness at significantly absorb vibrations and disturbances.





SOTO MODULAR Products Line - Vibration Absorption Systems

SDTDMD Products Line	Device Capabilities	Device Dimensions (mm)	Load Capacity (kg)
SOTOMOSOO	Vibration Absorption System	120x110x75 single isolator	Single Isolator can carry load up to 300kg.
SOTOMOSOO	Vibration Absorption System	120x110x75 single isolator	Single Isolator can carry load up to 500kg.
SOTOMO	Vibration Absorption System Stationed on multiple independent isolators	120x110x75 single isolator	Total load capacity for 4 isolators = SOTOMO X 4

SOTO MODULAR customized solutions available